

## CERTIFICATE OF ANALYSIS

# Certificate of Analysis

**Product name:** Modified Alsever's Solution

**Product code(s):** 04-182

**Batch Number(s):** 32180118

**Date Manufactured:** 23 Jan 18

**Expiry Date:** 22 Jan 19

**Appearance:** Clear Liquid

**Quality Control:**

	Specification	This batch result
pH	6.80-7.30	6.80
Conductivity( $\mu\text{s}/\text{cm}$ )	11500-13000	11830
Osmolality(mOsm/kg)	335-360	353

**Preservatives used:** Chloramphenicol and Neomycin Sulphate

**Sterility:** Terminally filtered with 0.2 $\mu\text{m}$  filter in a grade B clean room and tested for sterility.

**Storage:** 2-8oC

**Country of Origin:** United Kingdom

**Health and Safety:**

**Precautionary Notes**

Please download the Safety Data Sheet at <http://sourcebioscience.com>

**Signed for Quality Control:**



**Date:**

05 Feb 18

**Signed for Certificate Check:**



**Date:**

05 Feb 18