

CERTIFICATE OF ANALYSIS

Certificate of Analysis

Product name: Modified Alsever's Solution

Product code(s): 04-182, 04-182HE

Batch Number(s): 32180119

Date Manufactured: 29 Jan 18

Expiry Date: 28 Jan 19

Appearance: Clear Liquid

Quality Control:

	Specification	This batch result
pH	6.80-7.30	6.93
Conductivity($\mu\text{s}/\text{cm}$)	11500-13000	12650
Osmolality(mOsm/kg)	335-360	356

Preservatives used: Chloramphenicol and Neomycin Sulphate

Sterility: Terminally filtered with 0.2 μm filter in a grade B clean room and tested for sterility.

Storage: 2-8oC

Country of Origin: United Kingdom

Health and Safety:

Precautionary Notes

Please download the Safety Data Sheet at <http://sourcebioscience.com>

Signed for Quality Control: 

Date: 12 Feb 18

Signed for Certificate Check: 

Date: 14 Feb 18