

CERTIFICATE OF ANALYSIS

Certificate of Analysis

Product name: Red Cell Support Solution (RCSS)
Product code(s): 04-135
Batch Number(s): 35180110
Date Manufactured: 31 Jan 18
Expiry Date: 30 Jan 20
Appearance: Clear liquid
Quality Control:

	Specification	This batch result
pH	6.80-7.20	6.91
Conductivity(µs/cm)	3100-3300	3155
Osmolality(mOsm/kg)	300-315	311

Preservatives used: Antibiotics: Sulfamethoxazole and Trimethoprim

Sterility: Terminally filtered with 0.2µm filter in a grade B clean room and tested for sterility

Storage: 2-8oC

Country of Origin: United Kingdom

Health and Safety:

Precautionary Notes

Wear appropriate Personal Protective Equipment (PPE) when handling this reagent.
Please download the Safety Data Sheet at <http://sourcebioscience.com>

Signed for Quality Control:  **Date:** 12 Feb 18

Signed for Certificate Check:  **Date:** 21 Feb 18