

CERTIFICATE OF ANALYSIS

Certificate of Analysis

Product name: Heart Valve Cryopreservation Medium
Product code(s): 04-311
Batch Number(s): 71180208
Date Manufactured: 26 Feb 18
Expiry Date: 31 Oct 18
Appearance: Red translucent liquid
Quality Control:

	Specification	This batch result
pH	7.00-7.50	7.17
Conductivity($\mu\text{s}/\text{cm}$)	N/A	N/A
Osmolality(mOsm/kg)	N/A	N/A

Preservatives used: None

Sterility: Terminally filtered with 0.2 μm filter in grade B clean room and tested for sterility

Storage: Ambient

Country of Origin: United Kingdom

Health and Safety:

Precautionary Notes

Please download the Safety Data Sheet at <http://sourcebioscience.com>

Signed for Quality Control: 

Date: 12 Mar 18

Signed for Certificate Check: 

Date: 12 Mar 18