

CERTIFICATE OF ANALYSIS

Certificate of Analysis

Product name: Reagent Grade Water
Product code(s): 04-035
Batch Number(s): 18180385
Date Manufactured: 13 Mar 18
Expiry Date: 12 Mar 19
Appearance: Clear Liquid

Quality Control:

	Specification	This batch result
pH	N/A	N/A
Conductivity($\mu\text{s}/\text{cm}$)	<10	3.78
Osmolality(mOsm/kg)	N/A	N/A

e.g. Meets the requirements set out in 'Guidelines for the Blood Transfusion Services in the United Kingdom' 7th Edition; as set out in chapter 12, 12.2.4, Miscellaneous reagents.

Preservatives used: N/A

Sterility: Product is filtered through a 0.2 μm filter into gamma irradiated bags

Storage: Room temperature

Country of Origin: United Kingdom

Health and Safety:

Precautionary Notes

Please download the Safety Data Sheet at <http://sourcebioscience.com>

Signed for Quality Control: 

Date: 15 MAR 18

Signed for Certificate Check: 

Date: 15 March 18