

CERTIFICATE OF ANALYSIS

Certificate of Analysis

Product name: Tris Buffered Isotonic Saline (TBIS)

Product code(s): 04-142

Batch Number(s): 37190117

Date Manufactured: 08 Jan 19

Expiry Date: 30 Jun 20

Appearance: Clear Liquid

Quality Control:

	Specification	This batch result
pH	7.3 - 7.6	7.31
Conductivity($\mu\text{s}/\text{cm}$)	14,000 - 16,500	15096
Osmolality(mOsm/kg)	285 - 300	294

Manufactured to: SPN437/1 dated 11/05/12.

Preservatives used: None.

Sterility: Terminally filtered with a 0.2 μm filter in a grade B clean room and tested for sterility.e.g. Product is filtered through a 0.2 μm filter into gamma irradiated bags

Storage: Room temperature

Country of Origin: United Kingdom

Health and Safety:

Precautionary Notes

Please download the Safety Data Sheet at <http://sourcebioscience.com>

Signed for Quality Control:



Date: 14 Jan 19

Signed for Certificate Check:



Date: 14 JAN 19