

CERTIFICATE OF ANALYSIS

Certificate of Analysis

Product name: Modified Alsevers Solution
Product code(s): 05 - 182
Batch Number(s): 99190105
Date Manufactured: 14 Jan 19
Expiry Date: 12 Jan 20
Appearance: Clear Liquid

Quality Control:

	Specification	This batch result
pH	6.80 - 7.30	7.17
Conductivity($\mu\text{s}/\text{cm}$)	11500 - 13000	11596
Osmolality(mOsm/kg)	335 - 360	348

Preservatives used: Chloramphenicol and Neomycin sulphate

Sterility: Product is filtered through a 0.2 μm filter into gamma irradiated bags


Storage: 2 - 8°C

Country of Origin: United Kingdom

Health and Safety:

Precautionary Notes

Please download the Safety Data Sheet at <http://sourcebioscience.com>

Signed for Quality Control: 

Date: 21 JAN 19

Signed for Certificate Check: 

Date: 21 Jan 19