

CERTIFICATE OF ANALYSIS

Certificate of Analysis

Product name: Cambridge Antibiotic Solution
Product code(s): 04 - 301
Batch Number(s): 70190115
Date Manufactured: 14 Jan 19
Expiry Date: 27 Jul 19
Appearance: Red Suspension

Quality Control:

	Specification	This batch result
pH	7.1 - 7.5	7.26
Conductivity($\mu\text{s}/\text{cm}$)	N/A	N/A
Osmolality(mOsm/kg)	280 - 300	292

Manufacture Specification:
SPN310/2, effective 29/03/17

Preservatives used: Gentamicin sulphate, Imipenem, Polymyxin B, Vancomycin, Nystatin

Sterility: Product is filtered through a 0.2 μm filter into gamma irradiated bottles

Storage: <-40°C

Country of Origin: United Kingdom

Health and Safety:

Precautionary Notes

Please download the Safety Data Sheet at <http://sourcebioscience.com>

Signed for Quality Control: 

Date: 21 JAN 19

Signed for Certificate Check: 

Date: 21 Jan 19