

## CERTIFICATE OF ANALYSIS

# Certificate of Analysis

**Product name:** Cambridge Antibiotic Solution  
**Product code(s):** 04-301  
**Batch Number(s):** 70190216  
**Date Manufactured:** 04 Feb 19  
**Expiry Date:** 08 Feb 20  
**Appearance:** Red Suspension

**Quality Control:**

	Specification	This batch result
pH	7.1-7.5	7.10
Conductivity( $\mu\text{s}/\text{cm}$ )	N/A	N/A
Osmolality(mOsm/kg)	280-300	293

NHS BT Spec  
SPN310/2, effective 29/03/17.

**Preservatives used:** Gentamycin Sulphate, Imipenem, Polymixin B, Vancomycin, Nystatin

**Sterility:** Product is filtered through a 0.2 $\mu\text{m}$  filter prior to addition of Nystatin in a clean room and tested for micro sterility.

**Storage:** <-40 ° C

**Country of Origin:** United Kingdom

**Health and Safety:**

**Precautionary Notes**

Please download the Safety Data Sheet at <http://sourcebioscience.com>

**Signed for Quality Control:**



**Date:** 27 Feb 19

**Signed for Certificate Check:**



**Date:** 27 Feb 19