

CERTIFICATE OF ANALYSIS

Certificate of Analysis

Product name: Cambridge Antibiotic Solution
Product code(s): 04-301
Batch Number(s): 70190217
Date Manufactured: 18 Feb 19
Expiry Date: 31 Jan 20
Appearance: Red Suspension

Quality Control:

	Specification	This batch result
pH	7.1-7.5	7.3
Conductivity($\mu\text{s}/\text{cm}$)	N/A	N/A
Osmolality(mOsm/kg)	280-300	292

SPN310/2, effective 29/03/17.

Preservatives used: Gentamycin Sulphate, Imipenem, Polymixin B, Vancomycin, Nystatin

Sterility: Product is filtered through a 0.2 μm filter prior to addition of Nystatin in a clean room and tested for micro sterility.

Storage: <-40 ° C

Country of Origin: United Kingdom

Health and Safety:

Precautionary Notes

Please download the Safety Data Sheet at <http://sourcebioscience.com>

Signed for Quality Control:



Date: 28 Feb 19

Signed for Certificate Check:



Date: 28 FEB 19