

CERTIFICATE OF ANALYSIS

Certificate of Analysis

Product name: Modified Alsevers Solution
Product code(s): 05-182
Batch Number(s): 99190306
Date Manufactured: 11 Mar 19
Expiry Date: 16 Jan 20
Appearance: Clear Liquid
Quality Control:

	Specification	This batch result
pH	6.80-7.30	7.02
Conductivity($\mu\text{s}/\text{cm}$)	11500-13000	12297
Osmolality(mOsm/kg)	335-360	341

Preservatives used: Chloramphenicol, Neomycin sulphate and Levofloxacin

Sterility: Product is filtered through a 0.2 μm filter into gamma irradiated bags

Storage: 2 - 8°C

Country of Origin: United Kingdom

Health and Safety:

Precautionary Notes

Please download the Safety Data Sheet at <http://sourcebioscience.com>

Signed for Quality Control:



Date: 20 MAR 19

Signed for Certificate Check:



Date: 20 MAR 19