

## CERTIFICATE OF ANALYSIS

# Certificate of Analysis

**Product name:** Hypotonic Buffer Base

**Product code(s):** 04-333

**Batch Number(s):** 76190304

**Date Manufactured:** 20 Mar 19

**Expiry Date:** 16 Jan 20

**Appearance:** Clear Fluid

### Quality Control:

	Specification	This batch result
pH	7.80-8.20	8.09
Conductivity( $\mu\text{s}/\text{cm}$ )	N/A	N/A
Osmolality(mOsm/kg)	N/A	N/A

Manufactured to SPN320/2 Dated 29/03/17

**Preservatives used:** N/A

**Sterility:** Product is filtered through a 0.2 $\mu\text{m}$  filter and tested for micro sterility

**Storage:** 2-25°C

**Country of Origin:** United Kingdom

### Health and Safety:

#### Precautionary Notes

Please download the Safety Data Sheet at <http://sourcebioscience.com>

**Signed for Quality Control:** 

**Date:** 26 MAR 19

**Signed for Certificate Check:** 

**Date:** 26 Mar 19