

## CERTIFICATE OF ANALYSIS

# Certificate of Analysis

**Product name:** Fetal Cell Detection (FCD) Kit for Guest Medical  
**Product code(s):** 05-050 GL  
**Batch Number(s):** 05190352A/B  
**Date Manufactured:** 19 Mar 19  
**Expiry Date:** 23 Jan 20  
**Appearance:** Eluting: Black/Brown opaque liquid  
Counterstain: Pink/Red Translucent liquid

### Quality Control:

	Specification	This batch result
<b>pH</b>	Eluting: 1.80 - 2.10	Eluting: 2.07
<b>Conductivity(<math>\mu</math>s/cm)</b>	N/A	N/A
<b>Osmolality(mOsm/kg)</b>	N/A	N/A

This product has been assessed and Passed using the Kleihaur testing of negative and positive control samples. (Demonstration of fetal hemoglobin in erythrocytes of a blood smear). KLEIHAUR E, BRAUN H, BETKE K Klin Wochenschr 1957 (35):637-638

**Preservatives used:** Eluting: Alcohol based  
Counterstain: Sodium Azide (0.10% w/v)

**Sterility:** N/A

**Storage:** Ambient (Room Temperature)

**Country of Origin:** United Kingdom

### Health and Safety:

#### Precautionary Notes

Please download the Safety Data Sheet at <http://sourcebioscience.com>

**Signed for Quality Control:** 

**Date:** 03 Apr 19

**Signed for Certificate Check:** 

**Date:** 15 Apr 19