

CERTIFICATE OF ANALYSIS

Certificate of Analysis

Product name: Royal Brompton Heart Valve Antibiotic Solution (RBAS)
Product code(s): 04-371
Batch Number(s): 85190613
Date Manufactured: 17 Jun 19
Expiry Date: 01 Feb 20
Appearance: Red liquid

Quality Control:

	Specification	This batch result
pH	7.20 - 7.40	7.27
Conductivity($\mu\text{s}/\text{cm}$)	1300 - 1500	1453
Osmolality(mOsm/kg)	285 - 295	291

Preservatives used: N/A

Sterility: Terminally filtered with 0.2 μm filter in a grade B clean room and tested for sterility.

Storage: Store at -20°C

Country of Origin: United Kingdom

Health and Safety:

Precautionary Notes

Please download the Safety Data Sheet at <http://sourcebioscience.com>

Signed for Quality Control: 

Date: 25 Jun 19

Signed for Certificate Check: 

Date: 25 Jun 19