

## CERTIFICATE OF ANALYSIS

# Certificate of Analysis

**Product name:** Frozen Cell Recovery Solution

**Product code(s):** 04-161BAG, 04-162BAG

**Batch Number(s):** 41190622

**Date Manufactured:** 24 Jun 19

**Expiry Date:** 09 Aug 20

**Appearance:** Clear Liquid

**Quality Control:**

	Specification	This batch result
pH	6.40 - 6.80	6.67
Conductivity( $\mu\text{s}/\text{cm}$ )	44000 - 52000	46670
Osmolality(mOsm/kg)	1080 - 1120	1084

Meets the requirements set out in 'Guidelines for the Blood Transfusion Services in the United Kingdom' 7th Edition; as set out in chapter 12, 12.2.4, Miscellaneous reagents.

**Preservatives used:** EDTA (Anti-fungal)

**Sterility:** Product is filtered through a 0.2 $\mu\text{m}$  filter into gamma irradiated bags

**Storage:** 2-8°C

**Country of Origin:** United Kingdom

**Health and Safety:**

**Precautionary Notes**

Please download the Safety Data Sheet at <http://sourcebioscience.com>

**Signed for Quality Control:** 

**Date:** 11 JUL 19

**Signed for Certificate Check:** 

**Date:** 11 Jul 19