

## CERTIFICATE OF ANALYSIS

# Certificate of Analysis

**Product name:** Red Cell Support Solution  
**Product code(s):** 04 - 131  
**Batch Number(s):** 35190714  
**Date Manufactured:** 03 Jul 19  
**Expiry Date:** 21 Jul 20  
**Appearance:** Clear Liquid

**Quality Control:**

	Specification	This batch result
pH	6.80 - 7.20	7.10
Conductivity( $\mu\text{s}/\text{cm}$ )	3100 - 3300	3227
Osmolality(mOsm/kg)	300 - 315	305

**Preservatives used:** Antibiotics: Sulfamethoxazole and Trimethoprim

**Sterility:** Product is filtered through a 0.2 $\mu\text{m}$  filter into gamma irradiated bags

**Storage:** Store at 2 - 8°C

**Country of Origin:** United Kingdom

**Health and Safety:**

**Precautionary Notes**

Please download the Safety Data Sheet at <http://sourcebioscience.com>

**Signed for Quality Control:** 

**Date:** 11 JUL 19

**Signed for Certificate Check:** 

**Date:** 11 Jul 19