

CERTIFICATE OF ANALYSIS

Certificate of Analysis

Product name: Skin Cryopreservation Medium

Product code(s): 04 - 336

Batch Number(s): 73190810

Date Manufactured: 28 Aug 19

Expiry Date: 01 Sep 20

Appearance: Bright Red Liquid

Quality Control:

	Specification	This batch result
pH	7.10 - 7.50	7.35
Conductivity($\mu\text{s}/\text{cm}$)	N/A	N/A
Osmolality(mOsm/kg)	N/A	N/A

Manufactured to SPN 311 2

Preservatives used: N/A

Sterility: Product is filtered through a 0.2 μm filter into gamma irradiated bottles

Storage: 2 - 25 °C

Country of Origin: United Kingdom

Health and Safety:

Precautionary Notes

Please download the Safety Data Sheet at <http://sourcebioscience.com>

Signed for Quality Control: 

Date: 17 Sep 19

Signed for Certificate Check: 

Date: 17 Sep 19