

CERTIFICATE OF ANALYSIS

Certificate of Analysis

Product name: Fixol (Fixative Reagent for Fetal Cell Detection Kit)
Product code(s): 05-051
Batch Number(s): 51190958
Date Manufactured: 09 Sep 19
Expiry Date: 06 May 21
Appearance: Clear transparent liquid

Quality Control:

	Specification	This batch result
pH	8.50-9.50	8.74
Conductivity($\mu\text{s}/\text{cm}$)	1750-2150	1882
Osmolality(mOsm/kg)	N/A	N/A

This product has been assessed and Passed using the Kleihaur testing of negative and positive control samples. (Demonstration of fetal hemoglobin in erythrocytes of a blood smear). KLEIHAUR E, BRAUN H, BETKE K Klin Wochenschr 1957 (35):637-638

Preservatives used: Methanol based product

Sterility: N/A

Storage: Room temperature

Country of Origin: United Kingdom

Health and Safety:

Precautionary Notes

Please download the Safety Data Sheet at <http://sourcebioscience.com>

Signed for Quality Control: 

Date: 14 Oct 19

Signed for Certificate Check: 

Date: 14 Oct 19