

## CERTIFICATE OF ANALYSIS

# Certificate of Analysis

**Product name:** PBSS (Phosphate Buffered Serology Saline)  
**Product code(s):** 04 - 004  
**Batch Number(s):** 01191002  
**Date Manufactured:** 14 OCT 2019  
**Expiry Date:** 13 OCT 2021  
**Appearance:** Clear Liquid  
**Quality Control:**

|   | Specification | This batch result |
|---|---------------|-------------------|
| pH                                      | 6.95-7.05     | 6.95              |
| Conductivity( $\mu\text{s}/\text{cm}$ ) | N/A           | N/a               |
| Osmolality(mOsm/kg)                     | 285-295       | 287               |

Meets the requirements set out in 'Guidelines for the Blood Transfusion Services in the United Kingdom' 7th Edition; as set out in chapter 12, 12.2.4, Miscellaneous reagents.

**Preservatives used:** Sodium azide (0.01% [w/v])

**Sterility:** Product is filled in Gamma Irradiated bags

**Storage:** Room temperature


**Country of Origin:** United Kingdom

### Health and Safety:

#### Precautionary Notes


Please download the Safety Data Sheet at <http://sourcebioscience.com>

**Signed for Quality Control:**



**Date:** 16 Oct 2019

**Signed for Certificate Check:**



**Date:** 16 Oct 2019