

CERTIFICATE OF ANALYSIS

Certificate of Analysis

Product name: Modified Alsevers Solution

Product code(s): 04 - 182

Batch Number(s): 32191035

Date Manufactured: 17 Oct 19

Expiry Date: 16 Oct 20

Appearance: Clear Liquid

Quality Control:

	Specification	This batch result
pH	6.80 - 7.30	7.07
Conductivity($\mu\text{s}/\text{cm}$)	11500 - 13000	12103
Osmolality(mOsm/kg)	N/A	N/A

Manufactured to SPN870/2.1 Dated 18/08/16

Preservatives used: Chloramphenicol and Neomycin sulphate

Sterility: Product is filtered through a 0.2 μm filter into gamma irradiated bags

Storage: 2 - 8°C

Country of Origin: United Kingdom

Health and Safety:

Precautionary Notes

Please download the Safety Data Sheet at <http://sourcebioscience.com>

Signed for Quality Control: 

Date: 24 Oct 19

Signed for Certificate Check: 

Date: 24 Oct 2019