

## CERTIFICATE OF ANALYSIS

# Certificate of Analysis

**Product name:** Modified Alsevers Solution  
**Product code(s):** 05 - 182  
**Batch Number(s):** 99191008  
**Date Manufactured:** 09 Oct 19  
**Expiry Date:** 08 Oct 20  
**Appearance:** Clear Liquid  
**Quality Control:**

	Specification	This batch result
pH	6.80 - 7.30	6.90
Conductivity( $\mu\text{s}/\text{cm}$ )	11,500 - 13,000	11,653
Osmolality(mOsm/kg)	335 - 360	350

**Preservatives used:** Chloramphenicol, Neomycin Sulphate and Levofloxacin

**Sterility:** Product is filtered through a 0.2 $\mu\text{m}$  filter into gamma irradiated bags

**Storage:** Room temperature 2 - 8 °C

**Country of Origin:** United Kingdom

### Health and Safety:

#### Precautionary Notes

Please download the Safety Data Sheet at <http://sourcebioscience.com>

**Signed for Quality Control:** 

**Date:** 17 OCT 19

**Signed for Certificate Check:** 

**Date:** 17 Oct 19