

CERTIFICATE OF ANALYSIS

Certificate of Analysis

Product name: Modified Alsevers Solution

Product code(s): 04 - 182

Batch Number(s): 32191036

Date Manufactured: 22 Oct 19

Expiry Date: 21 Oct 20

Appearance: Clear Liquid

Quality Control:

	Specification	This batch result
pH	6.80 - 7.30	6.90
Conductivity($\mu\text{s}/\text{cm}$)	11500 - 13000	12203
Osmolality(mOsm/kg)	N/A	N/A

Manufactured to SPN 870/2.1 Dated 18.08.16

→ effective
6 Nov 18

Preservatives used: Chloramphenicol and Neomycin Sulphate

Sterility: Product is filtered through a 0.2 μm filter into gamma irradiated bags

Storage: 2 - 8°C

Country of Origin: United Kingdom

Health and Safety:

Precautionary Notes

Please download the Safety Data Sheet at <http://sourcebioscience.com>

Signed for Quality Control: *S.J. Taylor*

Date: 28 OCT 19

Signed for Certificate Check: *Waller*

Date: 28 OCT 2019