

CERTIFICATE OF ANALYSIS

Certificate of Analysis

Product name: Sodium Azide Free Phosphate Buffered Saline (pH7.3)
Product code(s): 04 - 047
Batch Number(s): 25191101
Date Manufactured: 11 Nov 19
Expiry Date: 12 May 21
Appearance: Clear Liquid
Quality Control:

	Specification	This batch result
pH	7.25 - 7.35	7.28
Conductivity($\mu\text{s}/\text{cm}$)	N/A	N/A
Osmolality(mOsm/kg)	285 - 295	285

Preservatives used: N/A

Sterility: Product is filtered through a 0.2 μm filter into gamma irradiated bottles


Storage: Room temperature

Country of Origin: United Kingdom

Health and Safety:

Precautionary Notes

Please download the Safety Data Sheet at <http://sourcebioscience.com>

Signed for Quality Control: 

Date: 15 NOV 19

Signed for Certificate Check: 

Date: 15 NOV 19