

CERTIFICATE OF ANALYSIS

Certificate of Analysis

Product name: Anti-Tn (Salvia Sclarea)
Product code(s): Component of 03-061
Batch Number(s): 49180310
Date Manufactured: 12 Mar 18
Expiry Date: 11 Mar 20
Appearance: Clear Yellow liquid

Quality Control:

	Specification	This batch result
pH	6.40 - 7.00	6.72
Conductivity($\mu\text{s}/\text{cm}$)	N/A	N/A
Osmolality(mOsm/kg)	285-295	285

Preservatives used: Sodium azide (0.1% [w/v])

Sterility: N/A

Storage: 2 - 8 °C

Country of Origin: United Kingdom

Health and Safety:

Precautionary Notes

Sodium azide is poisonous and is an irritant. It can be fatal if swallowed and may cause cardiac disturbances. Although present in these reagents in small quantities, appropriate health and safety interventions are recommended. Appropriate Personal Protective Equipment should be worn at all times whilst carrying out tests. Accidental exposure to the skin should be treated by washing with mild detergent and water for 15-20 minutes. Eyes should be washed using water or normal saline for the same duration.

Signed for Quality Control:



Date: 26 Mar 18

Signed for Certificate Check:



Date: 27 Mar 18