

## CERTIFICATE OF ANALYSIS

# Certificate of Analysis

**Product name:** PBS Azide free pH 7.00  
**Product code(s):** 04-101  
**Batch Number(s):** 06180401  
**Date Manufactured:** 03 Apr 18  
**Expiry Date:** 28 Aug 19  
**Appearance:** Clear Liquid

### Quality Control:

	Specification	This batch result
<b>pH</b>	6.95 - 7.05	6.98
<b>Conductivity(<math>\mu\text{s}/\text{cm}</math>)</b>	N/A	N/A
<b>Osmolality(mOsm/kg)</b>	285-295	291

**Preservatives used:** N/A

**Sterility:** Terminally filtered with 0.2 $\mu\text{m}$  filter in a grade B clean room and tested for sterility.

**Storage:** Ambient

**Country of Origin:** United Kingdom

### Health and Safety:

#### Precautionary Notes

Please download the Safety Data Sheet at <http://sourcebioscience.com>

**Signed for Quality Control:** 

**Date:** 16 Apr 18

**Signed for Certificate Check:** 

**Date:** 16 April 18