

CERTIFICATE OF ANALYSIS

Certificate of Analysis

Product name: Reagent Grade Water
Product code(s): 04-035
Batch Number(s): 18180487
Date Manufactured: 17 Apr 2018
Expiry Date: 16 Apr 2019
Appearance: Clear Liquid

Quality Control:

	Specification	This batch result
pH	N/A	N/A
Conductivity($\mu\text{s}/\text{cm}$)	<10	5.96
Osmolality(mOsm/kg)	N/A	N/A

Preservatives used: None

Sterility: No chemicals used in this preparation. Water is run through the reverse osmosis unit and collected without it being stored. It is collected straight into the appropriate bags or bottles with as little exposure to the atmosphere as possible


Storage: Room temperature

Country of Origin: United Kingdom

Health and Safety:

Precautionary Notes

Please download the Safety Data Sheet at <http://sourcebioscience.com>

Signed for Quality Control: 

Date: 18 APR 18.

Signed for Certificate Check: 

Date: 18 APR 18