

CERTIFICATE OF ANALYSIS

Certificate of Analysis

Product name: Isotonic Saline
Product code(s): 04-025 and 04-026
Batch Number(s): 66180401
Date Manufactured: 03 Apr 18
Expiry Date: 15 Oct 19
Appearance: Clear liquid
Quality Control:

	Specification	This batch result
pH	N/A	N/A
Conductivity($\mu\text{s}/\text{cm}$)	N/A	N/A
Osmolality(mOsm/kg)	285-295	292

Preservatives used: None

Sterility: Terminally filtered with 0.2 μm filter in grade B clean room and tested for sterility

Storage: Ambient

Country of Origin: United Kingdom

Health and Safety:

Precautionary Notes

Please download the Safety Data Sheet at <http://sourcebioscience.com>

Signed for Quality Control: 

Date: 18 Apr 18

Signed for Certificate Check: 

Date: 18 Apr 18