

CERTIFICATE OF ANALYSIS

Certificate of Analysis

Product name: Frozen Cell Recovery Solution (FCRS)

Product code(s): 04-161BAG, 04-162BAG

Batch Number(s): 41180111

Date Manufactured: 24 Jan 18

Expiry Date: 23 Jan 20

Appearance: Clear Liquid

Quality Control:

	Specification	This batch result
pH	6.40-6.80	6.56
Conductivity($\mu\text{s}/\text{cm}$)	44,000-52,000	47386
Osmolality(mOsm/kg)	1080-1120	1118

Preservatives used: EDTA (anti-fungal)

Sterility: Terminally filtered with 0.2 μm filter in a grade B clean room and tested for sterility.

Storage: 2-8oC

Country of Origin: United Kingdom

Health and Safety:

Precautionary Notes

This reagent contains a high concentration of sodium chloride but presents only a very minor risk to staff handling the solution

Signed for Quality Control:



Date: 05 Apr 18

Signed for Certificate Check:



Date: 05 Apr 18