

CERTIFICATE OF ANALYSIS

Certificate of Analysis

Product name: Polyagglutination Kit
Product code(s): 03-061
Batch Number(s): 79180412
Date Manufactured: 19 Apr 18
Expiry Date: 27 Feb 19
Appearance: Clear yellow/amber liquid(s)
Kit made up of four lectin extracts: 1. Arachis hypogea - 14180307, 2. Glycine soja - 15180307
3. Salvia Sclarea - 49180411, 4. Dolichos biflorus - 29180320

Quality Control:

	Specification	This batch result
pH	N/A	N/A
Conductivity($\mu\text{s}/\text{cm}$)	N/A	N/A
Osmolality(mOsm/kg)	N/A	N/A

No specific guidelines are laid down for a polyagglutination kit but the kit is expected to show a pattern of agglutination with affected human red cells which will classify them as being T, Tn, Tk or Cad. It may be possible to further determine Th and Tx affected red cells

Preservatives used: Sodium azide (0.1% [w/v])

Sterility: N/A

Storage: 2 - 8 °C

Country of Origin: United Kingdom

Health and Safety:

Precautionary Notes

Please download the Safety Data Sheet at <http://sourcebioscience.com>

Signed for Quality Control:



Date: 23 APR 18

Signed for Certificate Check:



Date: 23 Apr 18