

CERTIFICATE OF ANALYSIS

Certificate of Analysis

Product name: Detergent Buffer Base

Product code(s): 04-334

Batch Number(s): 77180501

Date Manufactured: 02 May 18

Expiry Date: 01 May 19

Appearance: Clear Colourless Liquid

Quality Control:

| | Specification | This batch result |
|-----------------------------------------|---------------|-------------------|
| pH | 7.80-8.20 | 7.94 |
| Conductivity($\mu\text{s}/\text{cm}$) | N/A | N/A |
| Osmolality(mOsm/kg) | N/A | N/A |

Meets the requirements set out in SPN321/2

Preservatives used: N/A

Sterility: Terminally filtered with 0.2 μm filter in a grade B clean room and tested for sterility

Storage: 2-25oC

Country of Origin: United Kingdom

Health and Safety:

Precautionary Notes

Please download the Safety Data Sheet at <http://sourcebioscience.com>

Signed for Quality Control:



Date:

18 MAY 18

Signed for Certificate Check:



Date:

18 MAY 18