

CERTIFICATE OF ANALYSIS

Certificate of Analysis

Product name: Modified Alsever's Solution

Product code(s): 04-182/04-182HE

Batch Number(s): 32180423

Date Manufactured: 30 Apr 18

Expiry Date: 29 Apr 19

Appearance: Clear Liquid

Quality Control:

	Specification	This batch result
pH	6.80-7.30	7.18
Conductivity($\mu\text{s}/\text{cm}$)	11500-13000	12160
Osmolality(mOsm/kg)	N/A	N/A

Preservatives used: Chloramphenicol and Neomycin Sulphate

Sterility: Terminally filtered with 0.2 μm filter in a grade B clean room and tested for sterility.

Storage: 2-8oC

Country of Origin: United Kingdom

Health and Safety:

Precautionary Notes

Please download the Safety Data Sheet at <http://sourcebioscience.com>

Signed for Quality Control:



Date: 21 MAY 18

Signed for Certificate Check:



Date: 21 MAY 18