

CERTIFICATE OF ANALYSIS

Certificate of Analysis

Product name: Cambridge Antibiotic Solution

Product code(s): 04-301

Batch Number(s): 70180510

Date Manufactured: 21 May 18

Expiry Date: 20 May 19

Appearance: Red Suspension

Quality Control:

	Specification	This batch result
pH	7.1 - 7.5	7.37
Conductivity($\mu\text{s}/\text{cm}$)	N/A	N/A
Osmolality(mOsm/kg)	280 - 300	290

Manufacture Specification:
SPN310/1, effective 29/03/17

Preservatives used: Gentamicin sulphate, Imipenem, Polymyxin B, Vancomycin, Nystatin

Sterility: Terminally filtered with 0.2 μm filter in a grade B clean room and tested for sterility.

Storage: <0°C

Country of Origin: United Kingdom

Health and Safety:

Precautionary Notes

Please download the Safety Data Sheet at <http://sourcebioscience.com>

Signed for Quality Control:



Date: 04 Jun 18

Signed for Certificate Check:



Date: 04 Jun 18