

CERTIFICATE OF ANALYSIS

Certificate of Analysis

Product name: Anti-A1 (Dolichos)5x Concentrate Blood Grouping Reagent

Product code(s): 03-042

Batch Number(s): 29180518

Date Manufactured: 04 Jun 18

Expiry Date: 23 May 20

Appearance: Yellow/Amber liquid

Quality Control:

	Specification	This batch result
pH	6.50-6.90	6.85
Conductivity($\mu\text{s}/\text{cm}$)	N/A	N/A
Osmolality(mOsm/kg)	500-750 (1:5*)	535

*Dilute 1 part Anti-A1 with 4 parts water

Avidity: Matches the minimum requirement to react within 2 minutes (23oC).

Preservatives used: Sodium azide (0.1% [w/v])

Sterility: N/A

Storage: 2 - 8 °C

Country of Origin: United Kingdom

Health and Safety:

Precautionary Notes

Please download the Safety Data Sheet at <http://sourcebioscience.com>

Signed for Quality Control: 

Date: 06 Jun 18

Signed for Certificate Check: 

Date: 06 Jun 18