

CERTIFICATE OF ANALYSIS

Certificate of Analysis

Product name: Reagent Grade Water

Product code(s): 04-035

Batch Number(s): 18180689

Date Manufactured: 06 Jun 2018

Expiry Date: 05 Jun 2019

Appearance: Clear Liquid

Quality Control:

| | Specification | This batch result |
|---|---------------|-------------------|
| pH | N/A | N/A |
| Conductivity($\mu\text{s}/\text{cm}$) | <10 | 2.85 |
| Osmolality(mOsm/kg) | N/A | N/A |

Preservatives used: None

Sterility: No chemicals used in this preparation. Water is run through the reverse osmosis unit and collected without it being stored. It is collected straight into the appropriate bags or bottles with as little exposure to the atmosphere as possible

Storage: Room temperature

Country of Origin: United Kingdom

Health and Safety:

Precautionary Notes

Please download the Safety Data Sheet at <http://sourcebioscience.com>

Signed for Quality Control: *D Crossto*

Date: 07-JUNE-18

Signed for Certificate Check: *on Q2*

Date: 07-JUNE 18