

CERTIFICATE OF ANALYSIS

Certificate of Analysis

Product name: Sterile 100% Glycerol
Product code(s): 04-361
Batch Number(s): 86180509
Date Manufactured: 24 May 18
Expiry Date: 15 June 20
Appearance: Translucent Viscous Liquid

Quality Control:

	Specification	This batch result
pH	5.5 - 8.0	7.27
Conductivity($\mu\text{s}/\text{cm}$)	N/A	N/A
Osmolality(mOsm/kg)	N/A	N/A

Complies with SPN434/1

Preservatives used: N/A

Sterility: Filled in Clean room and tested for sterility.

Storage: Room temperature

Country of Origin: United Kingdom

Health and Safety:

Precautionary Notes

Please download the Safety Data Sheet at <http://sourcebioscience.com>

Signed for Quality Control: 

Date: 11 Jun 18

Signed for Certificate Check: 

Date: 11 Jun 18