

## CERTIFICATE OF ANALYSIS

# Certificate of Analysis

**Product name:** LISS (Low Ionic Strength Saline)

**Product code(s):** 04-016

**Batch Number(s):** 03180622

**Date Manufactured:** 06 Jun 18

**Expiry Date:** 05 Jun 19

**Appearance:** Clear liquid

### Quality Control:

	Specification	This batch result
<b>pH</b>	6.50 - 6.90	6.81
<b>Conductivity(<math>\mu\text{s}/\text{cm}</math>)</b>	3500 - 3800	3605
<b>Osmolality(mOsm/kg)</b>	285 - 295	295

**Preservatives used:** Sodium azide (0.09% [w/v])

**Sterility:** N/A

**Storage:** Ambient

**Country of Origin:** United Kingdom

### Health and Safety:

#### *Precautionary Notes*

Wear appropriate Personal Protective Equipment (PPE) when handling this reagent.

Please download the Safety Data Sheet at <http://sourcebioscience.com>

**Signed for Quality Control:**



**Date:** 12 Jun 18

**Signed for Certificate Check:**



**Date:** 12 Jun 18