

CERTIFICATE OF ANALYSIS

Certificate of Analysis

Product name: PBSS AF (Phosphate Buffered Serology Saline, Azide Free) pH 7.30

Product code(s): 04-047

Batch Number(s): 25180508

Date Manufactured: 30 May 18

Expiry Date: 27 Aug 19

Appearance: Clear Liquid

Quality Control:

	Specification	This batch result
pH	7.25 - 7.35	7.33
Conductivity($\mu\text{s}/\text{cm}$)	N/A	N/A
Osmolality(mOsm/kg)	285 - 295	295

Preservatives used: N/A

Sterility: Terminally filtered with a 0.2 μm sterile filter in a grade B clean room and tested for sterility.

Storage: Room Temperature

Country of Origin: United Kingdom

Health and Safety:

Precautionary Notes

Please download the Safety Data Sheet at <http://sourcebioscience.com>

Signed for Quality Control: 

Date: 21 Jun 18

Signed for Certificate Check: 

Date: 21 Jun 18