

CERTIFICATE OF ANALYSIS

Certificate of Analysis

Product name: AmnioCryo Medium
Product code(s): 04-326
Batch Number(s): 72180511
Date Manufactured: 17 May 18
Expiry Date: 14 Mar 20
Appearance: Clear Red Solution

Quality Control:

	Specification	This batch result
pH	7.10 - 7.50	7.31
Conductivity($\mu\text{s}/\text{cm}$)	N/A	N/A
Osmolality(mOsm/kg)	N/A	N/A

Manufactured to NHSBT Specification: SPN431, dated 06/07/17.

Preservatives used: None

Sterility: Terminally filtered with 0.2 μm filter in a grade B clean room and tested for sterility.

Storage: Ambient

Country of Origin: United Kingdom

Health and Safety:

Precautionary Notes

Please download the Safety Data Sheet at <http://sourcebioscience.com>

Signed for Quality Control:



Date: 09 Jul 18.

Signed for Certificate Check:



Date: 10 Jul 18