

CERTIFICATE OF ANALYSIS

Certificate of Analysis

Product name: Cryoprotectant Reagent

Product code(s): 04-213

Batch Number(s): 61180710

Date Manufactured: 19 Jul 18

Expiry Date: 04 Apr 20

Appearance: Clear Liquid

Quality Control:

	Specification	This batch result
pH	3.50 - 5.50	4.75
Conductivity($\mu\text{s}/\text{cm}$)	4200-4850	4527
Osmolality(mOsm/kg)	880-970	955

Preservatives used: None

Sterility: Terminally filtered with 0.2 μm filter in a grade B clean room and tested for sterility.

Storage: Ambient

Country of Origin: United Kingdom

Health and Safety:

Precautionary Notes

Please download the Safety Data Sheet at <http://sourcebioscience.com>

Signed for Quality Control: 

Date: 26 Jul 18

Signed for Certificate Check: 

Date: 26 Jul 18