

## CERTIFICATE OF ANALYSIS

# Certificate of Analysis

**Product name:** Reagent Grade Water

**Product code(s):** 04-035

**Batch Number(s):** 18180892

**Date Manufactured:** 02 Aug 2018

**Expiry Date:** 01 Aug 2019

**Appearance:** Clear liquid

**Quality Control:**

	Specification	This batch result
pH	N/A	N/A
Conductivity( $\mu\text{s}/\text{cm}$ )	0 - 10	0.085
Osmolality(mOsm/kg)	N/A	N/A

**Preservatives used:** No preservatives used

**Sterility:** Water is UV treated and processed into gamma irradiated bags. Although exposure to the environment is limited, we do not guarantee the sterility of the water.

**Storage:** Store at Room Temperature

**Country of Origin:** United Kingdom

**Health and Safety:**

**Precautionary Notes**

There are no H&S implications from this product;

Please download the safety data sheet from [www.sourcebioscience.com](http://www.sourcebioscience.com)

**Signed for Quality Control:** 

**Date:** 06-Aug-18

**Signed for Certificate Check:**



**Date:** 06 Aug 18