

CERTIFICATE OF ANALYSIS

Certificate of Analysis

Product name: Modified Alsevers Solution

Product code(s): 04-182

Batch Number(s): 32180725

Date Manufactured: 30 Jul 18

Expiry Date: 29 Jul 19

Appearance: Clear Liquid

Quality Control:

	Specification	This batch result
pH	6.80 - 7.30	6.93
Conductivity($\mu\text{s}/\text{cm}$)	11,500 - 13,000	1173
Osmolality(mOsm/kg)	N/A	N/A

Manufacture Specification:
SPN870/1, effective 18/08/16.

Preservatives used: Chloramphenicol and Neomycin Sulphate

Sterility: Terminally filtered with 0.2 μm filter in a grade B clean room and tested for sterility.

Storage: 2-8°C

Country of Origin: United Kingdom

Health and Safety:

Precautionary Notes


Please download the Safety Data Sheet at <http://sourcebioscience.com>

Signed for Quality Control:



Date: 22 Aug 18

Signed for Certificate Check:



Date: 22 Aug 18