

CERTIFICATE OF ANALYSIS

Certificate of Analysis

Product name: LISS 2 (LISS Additive Reagent)
Product code(s): 04-070
Batch Number(s): 34180812
Date Manufactured: 21 Aug 18
Expiry Date: 04 Apr 20
Appearance: Clear Liquid

Quality Control:

	Specification	This batch result
pH	6.70-7.10	7.03
Conductivity($\mu\text{s}/\text{cm}$)	3000-3600	3153
Osmolality(mOsm/kg)	290-310	296

Meets the requirements set out in 'Guidelines for the Blood Transfusion Services in the United Kingdom' 7th Edition; as set out in chapter 12, 12.2.4, Miscellaneous reagents.

Preservatives used: Sodium azide (0.08% [w/v])

Sterility: Terminally filtered with 0.2 μm filter in a grade B clean room and tested for sterility.

Storage: 2-8oC

Country of Origin: United Kingdom

Health and Safety:

Precautionary Notes

Please download the Safety Data Sheet at <http://sourcebioscience.com>

Signed for Quality Control:



Date:

13 Sep 18

Signed for Certificate Check:



Date:

13 Sep 18