

CERTIFICATE OF ANALYSIS

Certificate of Analysis

Product name: Tris Buffered Isotonic Saline (TBIS)
Product code(s): 04-141 and 04-142
Batch Number(s): 37180611
Date Manufactured: 07 Jun 18
Expiry Date: 21 Mar 19
Appearance: e.g. Clear Liquid

Quality Control:

	Specification	This batch result
pH	7.3 - 7.6	7.43
Conductivity($\mu\text{s}/\text{cm}$)	14,000 - 16,500	15130
Osmolality(mOsm/kg)	285 - 300	295

Manufactured to: SPN437/1 dated 11/05/12.

Preservatives used: None.

Sterility: Terminally filtered with a 0.2 μm filter in a grade B clean room and tested for sterility.e.g. Product is filtered through a 0.2 μm filter into gamma irradiated bags

Storage: e.g. Room temperature

Country of Origin: United Kingdom

Health and Safety:

Precautionary Notes

Please download the Safety Data Sheet at <http://sourcebioscience.com>

Signed for Quality Control:



Date: 18 Sep 18.

Signed for Certificate Check:



Date: 18 Sep 18