

## CERTIFICATE OF ANALYSIS

# Certificate of Analysis

**Product name:** Fixol (Fixative reagent for Fetal Cell Detection Kit) for guest medical

**Product code(s):** 05-051

**Batch Number(s):** 51180942C

**Date Manufactured:** 13 Sep 18

**Expiry Date:** 22 Apr 20

**Appearance:** Clear transparent liquid

**Quality Control:**

	Specification	This batch result
pH	8.50 - 9.50	8.89
Conductivity( $\mu\text{s}/\text{cm}$ )	1750 - 2150	1770
Osmolality(mOsm/kg)	N/A	N/A

This product has been assessed and Passed using the Kleihaur testing of negative and positive control samples. (Demonstration of fetal hemoglobin in erythrocytes of a blood smear). KLEIHAUR E, BRAUN H, BETKE K Klin Wochenschr 1957 (35):637-638

**Preservatives used:** e.g. Sodium azide (0.01% w/v), Alcohol based

**Sterility:** N/A

**Storage:** Ambient (Room Temperature)

**Country of Origin:** United Kingdom

**Health and Safety:**

### Precautionary Notes

Please download the Safety Data Sheet at <http://sourcebioscience.com>

**Signed for Quality Control:**



**Date:** 15 Oct 18

**Signed for Certificate Check:**



**Date:** 15 Oct 18