

CERTIFICATE OF ANALYSIS

Certificate of Analysis

Product name: Heart Valve Cryopreservation Medium
Product code(s): 04-311
Batch Number(s): 71181010
Date Manufactured: 01 Oct 18
Expiry Date: 30 Oct 19
Appearance: Clear Red Solution

Quality Control:

	Specification	This batch result
pH	7.00 - 7.50	7.44
Conductivity($\mu\text{s}/\text{cm}$)	NA	NA
Osmolality(mOsm/kg)	NA	NA

NHSBT Specification
SPN323/2, Dated 24/09/18.

Preservatives used: None

Sterility: Product is filtered through a 0.2 μm filter into gamma irradiated bottles

Storage: 0 - 30°C

Country of Origin: United Kingdom

Health and Safety:

Precautionary Notes

Please download the Safety Data Sheet at <http://sourcebioscience.com>

Signed for Quality Control:



Date: 22 Oct 18

Signed for Certificate Check:



Date: 22 Oct 18