

CERTIFICATE OF ANALYSIS

Certificate of Analysis

Product name: LISS (Low Ionic Strength Saline)
Product code(s): 04-016, 04-123
Batch Number(s): 03181026
Date Manufactured: 22 Oct 18
Expiry Date: 13 Aug 19
Appearance: Clear Liquid
Quality Control:

	Specification	This batch result
pH	6.50-6.90	6.77
Conductivity($\mu\text{s}/\text{cm}$)	3500-3800	3628
Osmolality(mOsm/kg)	285-295	285

Preservatives used: Sodium azide (0.09% [w/v])

Sterility: N/A

Storage: Ambient

Country of Origin: United Kingdom

Health and Safety:

Precautionary Notes

Wear appropriate Personal Protective Equipment (PPE) when handling this reagent.
Please download the Safety Data Sheet at <http://sourcebioscience.com>

Signed for Quality Control:



Date: 29 Oct 18

Signed for Certificate Check:



Date: 30 Oct 18