

## CERTIFICATE OF ANALYSIS

# Certificate of Analysis

**Product name:** Detergent Buffer Base

**Product code(s):** 04-334

**Batch Number(s):** 77181003

**Date Manufactured:** 18 Oct 18

**Expiry Date:** 16 Jan 20

**Appearance:** Clear Liquid

### Quality Control:

	Specification	This batch result
pH	7.8-8.2	8.19
Conductivity( $\mu\text{s}/\text{cm}$ )	N/A	N/A
Osmolality(mOsm/kg)	N/A	N/A

NHS BT Spec

Meets the requirements set out in SPN321/2.

**Preservatives used:** N/A

**Sterility:** Terminally filtered with 0.2 $\mu\text{m}$  filter in a grade C clean room and tested for sterility.

**Storage:** 2-25 °C

**Country of Origin:** United Kingdom

### Health and Safety:

#### Precautionary Notes

Please download the Safety Data Sheet at <http://sourcebioscience.com>

**Signed for Quality Control:**



**Date:**

31 Oct 18

**Signed for Certificate Check:**



**Date:**

31 Oct 18