

CERTIFICATE OF ANALYSIS

Certificate of Analysis

Product name: Nuclease Buffer Base

Product code(s): 04-335

Batch Number(s): 78181004

Date Manufactured: 17 Oct 18

Expiry Date: 30 Jun 20

Appearance: Clear Liquid

Quality Control:

	Specification	This batch result
pH	7.3-7.7	7.68
Conductivity($\mu\text{s}/\text{cm}$)	N/A	N/A
Osmolality(mOsm/kg)	N/A	N/A

NHS BT Spec
SPN322/2 Refers effective 28/03/17

Preservatives used: N/A

Sterility: Product is filtered through a 0.2 μm filter in grade B clean room and tested for micro sterility.

Storage: 2-25 ° C

Country of Origin: United Kingdom

Health and Safety:

Precautionary Notes


Please download the Safety Data Sheet at <http://sourcebioscience.com>

Signed for Quality Control:



Date: 31 Oct 18

Signed for Certificate Check:



Date: 31 Oct 18