

## CERTIFICATE OF ANALYSIS

# Certificate of Analysis

**Product name:** Hydrochloric Acid 0.5M

**Product code(s):** 04-356

**Batch Number(s):** 84181012

**Date Manufactured:** 24 Oct 18

**Expiry Date:** 23 Oct 20

**Appearance:** Clear Liquid

**Quality Control:**

	Specification	This batch result
pH	N/A	N/A
Conductivity( $\mu$ s/cm)	N/A	N/A
Osmolality(mOsm/kg)	920-1020	996

NHS BT Spec  
Referred to SPN435/1

**Preservatives used:** N/A

**Sterility:** Product is filtered through a 0.2 $\mu$ m filter and tested for micro sterility.

**Storage:** Ambient

**Country of Origin:** United Kingdom

**Health and Safety:**

**Precautionary Notes**

Please download the Safety Data Sheet at <http://sourcebioscience.com>

**Signed for Quality Control:**



**Date:**

01 NOV 18

**Signed for Certificate Check:**



**Date:**

01 NOV 18