

CERTIFICATE OF ANALYSIS

Certificate of Analysis

Product name: AmnioCryo Medium
Product code(s): 04-326
Batch Number(s): 72181011
Date Manufactured: 25 Oct 18
Expiry Date: 13 Mar 20
Appearance: Clear Red Solution

Quality Control:

	Specification	This batch result
pH	7.10-7.50	7.11
Conductivity($\mu\text{s}/\text{cm}$)	N/A	N/A
Osmolality(mOsm/kg)	N/A	N/A

NHS BT Spec

Manufactured to SPN431/2, dated 06/07/17.

Preservatives used: N/A

Sterility: Product is filtered through a 0.2 μm filter and tested for micro sterility.

Storage: Ambient

Country of Origin: United Kingdom

Health and Safety:

Precautionary Notes

Please download the Safety Data Sheet at <http://sourcebioscience.com>

Signed for Quality Control:



Date:

06 Nov 18

Signed for Certificate Check:



Date:

06 Nov 18